



PRODUCT COMPARISON

PERFORMANCELINE

UFSK
INTERNATIONAL

SYS



ophthalmic surgery



IVOM treatments



implantology
maxillofacial surgery



hand surgery



500 XLE

Item No.: 1500 0006



500 XLS

Item No.: 1500 0005



chassis / dimensions

805 x 605 mm

1010 x 710 mm



material

zinc undercoated, powder-coated
welded steel frame

painting steel frame

casters

4 big double casters
4x 150 mm, 360°
braking, neutral
steering,
electrically conducting

4 casters
2x 100 mm double casters, 360°
braking
2x 150 mm, fixed
steering

braking system

manual

central brake

No, brake left / right

No, brake front casters

brake lever

left / right foot lever

steering

Yes (manual)

Yes

foot switch

removable foot lever

foot switch functions

on / off (POWER)
adjusting head, height

on / off (POWER)
adjusting back, seat., height
3-memo- / reset-position

top frame

zinc undercoated, powder-coated
welded steel frame

zinc undercoated,
painting steel frame

SES - surface extension system

No

No

bed surface / wide

600 mm

norm rails

Yes (optional)
Item No.: 3351 0047

Yes (optional)
Item No.: 3321 0000 / 3320 0000

upholstery

back 72 mm
seat 75 mm
leg 60 mm

back 100 mm
seat 100 mm
leg 80 mm

„Tundra skai“

headrest

foam PUR, seamless, with hollow



U shape, seams, small hollow



headpart / adjustments

electrical, manual

manual

pushing handles

removable, stainless steel with rubber

fixed, stainless steel



Distribuido por: Accutomex SA de CV

www.accutomex.com

ventas@accutomex.com

Tel. +52 (81) 8347-2520 / 01-800-800-0507



PRODUCT COMPARISON

PERFORMANCELINE

UFSK
INTERNATIONAL

SYS



ophthalmic surgery



IVOM treatments



implantology
maxillofacial surgery



hand surgery




500 XLE

Item No.: 1500 0006



500 XLS

Item No.: 1500 0005



Hand control functions	on / off (POWER) adjusting head, back, seat, height	on / off (POWER) adjusting back, seat, height
battery display LED		No
brake activ LED		No
steering guide activ LED		No
reset position		Yes
memo positionen		4
max. weight	250 kg	200 kg
battery	back up battery Item No.: 3047 0000	
mains operation	No	
surgery cycle	aprox. 100	
e-drive	No	
entry height	600 mm (standard) optional 520 mm with foot extension Item No.: 3351 0049	580 mm (standard)