



eyeForce



500 XLE



## Our quality features at a single glance

In addition to bringing you an optimized price-performance ratio, all our products come with special quality features designed to benefit staff, surgeons, and patients alike while making daily work life easier, especially in high patient-turnover surgery environments.



**Advantages  
for the surgeon**



**Advantages  
for the assistant**



**Advantages  
for the patient**



### Steering Guide

- Facilitates the directional pushing by only a single person while protecting that person's back muscles against strain



### Mobility

- Its lightweight construction reduces the treatment chair's weight
- 360-degree rotating double casters



### Memory Auto-Run

- Move chair into pre-programmed positions at the push of a button
- The auto-run functions moves the chair automatically into the desired position
- Ergonomic design to protect the back of staff and patients alike



### Electronic brake system

- Allows you to set the brake with the push of a button – for no more cumbersome brake levers in the staff's foot area

500 XLE comfort



500 XLS



500 bariatric



**Max. payload**

- Maximum patient weight admissible on the equipment



**Ergonomics**

- Optimized patient positioning
- Ergonomically adjustable 4-element upper frame
- All four elements of the support surface are separately adjustable via the manual keypad



**Pop-up joystick**

- Foot operation via joystick for sensitive chair adjustment



**Battery powered**

- Battery-operated; independent from the power supply; perfectly mobile
- Powerful, eco-friendly batteries
- User-friendly back-up system



**Premium upholstery**

- High-grade, easy-to-clean, and shape retaining
- Anti-microbially treated, hygienic, mildew-retardant
- Flame retardant and waterproof



**e-drive system**

- Electronic pushing aid, very easy to use



**Enormous legroom**

- Ensure optimum surgeon's access to the patient



**Headrest element**

- 3-dimensionally adjustable headrest element
- Manually or electronically adjustable



**Extension system SES**

- Low patient access height and very little need for repositioning thanks to automatic length compensation by the footrest and legrest elements
- The patient remains completely relaxed



**Foot keypad**

- Removable foot keypad allowing a variety of adjusting options



**Service-friendly**

- Easy maintenance and care