

SL7D

110V/220V Haag Strait Type 5 Magnifications Tonometer
adaptor included Yellow (Yellow or White LED illumination)


Häns Heiss

MICROSCOPE

Type	Galilean-Type
Configuration	Tower style
Magnification change	Five Position revolving Drum
Eyepieces	12.5X
Angle between eyepieces	13°
Total magnification Ratio	6X, 10X, 16X, 25X, 40X
Pupillary adjustment	52mm~78mm
Diopter adjustment	±6D
Field of view	40X (Ø5.5mm), 25X (Ø8.5mm), 16X (Ø13.5mm), 10X (Ø22mm), 6X (Ø34.7mm)

SLIT ILLUMINATION

Slit width	Continuously variable from 0 to 14mm (at 14mm, slit becomes a circle)
Slit length	Continuously variable from 1mm to 14mm
Aperture diameters	Ø14mm. Ø10mm. Ø5mm. Ø3mm, Ø2mm. Ø1mm, Ø0.2mm
Slit angle	0°-180°
Slit inclination 4 step	5°, 10°, 15°, 20°
Filters	thermal safety, UV, Red-free, Cobalt Blue
Lamp	6V/20W Halogen Lamp (OSRAM, Import from Germany)



Distribuido por: Accutomex SA de CV

www.accutomex.com

ventas@accutomex.com

Tel. +52 (81) 8347-2520 / 01-800-800-0507

SL7D
**110V/220V Haag Strait Type 5 Magnifications Tonometer
adaptor included Yellow (Yellow or White LED illumination)**
Häns Heiss

BASE

Longitudinal movement	90mm
Lateral movement	100mm
Fine Base movement	15mm
Vertical movement	30mm

Chin-Rest

Vertical movement	80mm
Fixation Target	Red LED

Power

Input voltage	220V/110V~±10%
Input frequency	50Hz/60Hz
Power Consumption	30VA (max)

Output Voltage

Light	6V
Fixation	3V

Dimension & Weight

Dimension	740mm x 450mm x 500mm
Gross weight	25Kg
Net weight	24Kg

Working Environment

Temperature	+5°C~+40°C
Relative humidity	≤80%
Air pressure	800hpa~1060hpa

Storing Environment

Temperature	-40°C~+55°C
Relative humidity	≤93%
Air pressure	700hpa~1060hpa

Transporting Environment

Temperature	-40°C~+55°C
Relative humidity	≤93%
Air pressure	700hpa~1060hpa

VIDEO SYSTEM B was upgraded by the Slit Lamps with the CCD Adaptor, Beam Splitter and CCD with high speed and resolution. Besides the functions of the Slit Lamps it also has the following functions: 1) 5.0 mega high-definition static image capture. 2) Easy operated patients' information management, diagnosis report can be edited.

ACCESSORIES Applanation Tonometer, Beam Splitter, DC Adaptor, DC, Video Camera Adaptor, Video Camera, BacLite Background Illumination Model, Teaching Tube, etc... Yellow Filter, Amber Filter, Polarized Filter and Measuring Eyepiece are available.

SL7D

110V/220V Haag Strait Type 5 Magnifications Tonometer
adaptor included Yellow (Yellow or White LED illumination)

Häns Heiss

