

surgiForce

The ultimate surgeon's chair for microsurgery.



MICRO
SURGERY



OPHTHALMIC
SURGERY

Order No. 26000000

ADVANTAGES



Armrests

- Allows personalized adjustment of the armrest elements under sterile conditions even while performing surgical procedures - simply via the foot keypad and the single-handed activation

Seating system

- Our unique dynamic seating system supports correct sitting postures while protecting your back
- The seat surface and seat depth can be individually adjusted to suit every surgeon's physique with the aid of a sliding carrier plate (located underneath the seat).
- Moving the seat allows individualized ergonomic adjustment of the surgical area and ensures that the patient is sufficiently close to the attending physician
- Seat height adjustment is also effected electronically, via the foot switches


Chassis

- The cast-aluminum lightweight chassis featuring double casters and a central brake provides maximum stability and at the same time unlimited mobility
- The materials and workmanship satisfy the most stringent hygienic requirements
- Battery operation for power supply-independent use
- Foot control element for seat height and armrest height adjustment; acoustic signal tone signaling the switch-over from armrest height to seat height adjustment
- Added stability thanks to a 5th wheel



surgiTrend

*Designed with preventive medicine
in mind.*



Order No. 22000010

ADVANTAGES



Armrests

- Rugged, 360-degree movable, and height-adjustable

Seat System

- Our unique dynamic seating system supports correct sitting postures while protecting your back

Chassis

- 4 double casters
- 5th wheel for added resistance to tipping
- Central brake acting on 4 wheels
- Power supply-independent battery box with integrated emergency stop button
- Foot keypad panel for moving the seat up & down



surgiLine

The modern surgeon's stool.



OPHTHALMIC
SURGERY

Order No. 22000008

ADVANTAGES



Armrests

- Rugged, 360-degree movable, and height-adjustable
- Designed for personalized adjustment to the surgeon

Seat System

- Seat system evenly distributes your weight as you make the slightest adjustment in position.
- Anatomically formed, ergonomic backrest reduces pressure on inter vertebral discs and provides ideal support for the lumbar spine.
- Seat system upholstery in standard color anthracite F6461197
- Outstanding sitting comfort

Chassis

- 4 double casters
- Central brake acting on 4 wheels
- Power supply-independent battery box with integrated emergency stop button
- Foot keypad for moving the seat up & down
- Extremely easy to move
- Space-saving



assistTrend

The ergonomic assistant's stool.



OPHTHALMIC
SURGERY



AESTHETIC
SURGERY



DENTAL
SURGERY



ENT
SURGERY



MICRO
SURGERY



ORAL &
MAXILLOFACIAL

Order No. 25000002

ADVANTAGES



- Our unique dynamic seating system supports correct sitting postures while protecting your back
- High-grade, comfortable and easy-care upholstery fabric to meet even the most stringent hygiene requirements
- Seat and backrest upholstery in standard color black
- Upholstery optionally available in a wide variety of colors
- Excellent movement characteristics thanks to double castors and 5-arm aluminum chassis
- Sterile, comfortable height adjustment via foot-triggered gas pressure spring
- Flexible height adjustment from 57 to 77 centimeters
- Reduces muscular tension in the neck and lumbar areas caused by work-related, prolonged sitting.
- Supports maximum loads up to 130 kilograms



Comfortable height adjustment
via foot triggering.



Optional: assistTrend with armrests.

No. 33000000



SPECIFICATIONS

Chassis		Max. Load	130 kg
- Five legged chassis	ø 585 mm		
- Wheels	ø 65 mm	Seat/backrest	
Lifting Column		- Seat/Back adjusting	-6°/+3°
- Height adjustment	57–77 cm	Upholstery	Stamskin top
- For body height of 1.75 m–1.95 m		Standard color	black
Weight	7,5 kg	Subject to technical modifications.	

assistTrend





OPHTHALMIC
SURGERY



AESTHETIC
SURGERY



DENTAL
SURGERY



ENT
SURGERY



MICRO
SURGERY




ORAL &
MAXILLOFACIAL

Order No. 24000015

ADVANTAGES



- Rugged safety chassis with 5-arm aluminum pedestal
- Double casters for extremely smooth movement
- Height adjustment via manually triggered gas pressure spring
- Flexible height adjustment from 54.5 to 73.5 centimeters
- Rotating and tilting anatomically correct saddle seat, comfortable and individually adjustable
- Maximum comfort level for your assistants, to prevent pain and fatigue while supporting the blood circulation
- High-quality, comfortable and easy-care upholstery to meet even the most stringent demands and hygiene requirements
- Seat upholstery in standard color anthracite (Color chart Tundra skai)
- Upholstery optionally available in a wide variety of colors
- Supports maximum loads up to 130 kilograms



Manual gas pressure spring
adjustment.



SPECIFICATIONS

Lifting Column

- Hand operated 190 mm
- Height adjustment 54.5-73.5 cm

Chassis

- Five legged chassis ø 585 mm
- Wheels ø 65 mm

Weight

Saddle seat 35 x 42.5 x 18 cm 7,5 kg


Seat

- Seat adjusting +2°/+16°

Upholstery

Tundra skai
Standard color anthracite

Subject to technical modifications.



assistPro

The perfect assistant's stool.



OPHTHALMIC
SURGERY



AESTHETIC
SURGERY



DENTAL
SURGERY



ENT
SURGERY



MICRO
SURGERY



ORAL &
MAXILLOFACIAL

Order No. 24000005 & 24000009

ADVANTAGES



- Ergonomically shaped seat and backrest
- The upholstery fabric is easy to clean, resistant to blood, perspiration and saliva, non-fade, and flame-retardant
- Excellent movement characteristics thanks to double castors and 5-arm aluminum chassis
- Sterile, comfortable height adjustment via **foot-triggered** gas pressure spring for the **assistPro F**. The height of **assistPro H** is adjusted **manually**, via the levers underneath the seat.
- Flexible height adjustment from 52 to 71.5 centimeters
- Reduces muscular tension in the neck and lumbar areas caused by work-related, prolonged sitting.
- Supports maximum loads up to 130 kilograms
- Upholstery in standard color black (Color chart Periline)


assistPro F

Order No. 24000005



assistPro H

Order No. 24000009



*Optional:
assistPro F with arm rests
Rugged armrests, 360-de-
gree movable and
height-adjustable
Order No.: 33510038*



*Optional:
Electrically conductive
upholstery in color black
Order No.: 30370223*

SPECIFICATIONS

Chassis

- Five legged chassis ø 585 mm
- Wheels ø 65 mm

Weight

- assistPro F 12 kg
- assistPro H 11 kg

Lifting Column

- Height adjustment 52–71.5 cm

Seat/backrest

- Seat adjusting +6°/-18°
- Back adjusting +4°/-9°
- Max. Load 130 kg
- Upholstery Periline
- Standard color black

Subject to technical modifications.

assistPro



easyPro

Versatile and space saving.



OPHTHALMIC
SURGERY



AESTHETIC
SURGERY



DENTAL
SURGERY



ENT
SURGERY



MICRO
SURGERY



ORAL &
MAXILLOFACIAL

Order No. 24000010, 24000006 & 24000013

ADVANTAGES



- Ergonomically shaped seat
- The upholstery fabric is easy to clean, resistant to blood, perspiration and saliva, non-fade, and flame-retardant
- Excellent movement characteristics thanks to double castors and 5-arm aluminum chassis
- Comfortable height adjustment via gas pressure spring with **foot-triggering** (sterile) for the **easyPro F**, with **manual triggering** for the **easyPro H**, and with **release-ring triggering** for the **easyPro R**
- Flexible height adjustment: easyPro F and easyPro H from 50 to 70 centimeters and easyPro R from 52 to 70 centimeters
- Reduces muscular tension in the neck and lumbar areas caused by work-related, prolonged sitting.
- Supports maximum loads up to 130 kilograms
- Upholstery in standard color black (Color chart Periline)



easyPro F
Order No. 24000010



easyPro H
Order No. 24000006



easyPro R
Order No. 24000013



*Optional:
Upholstery available
in a wide variety of
colors (Color chart
Periline page 69)*



*Optional:
Double rimmed castor,
65mm, with braking
function
Order No.: 33510040*

SPECIFICATIONS

Chassis		Seat	
- Five legged chassis	ø 585 mm	- Seat adjusting	-12°
- Wheels	ø 65 mm		
Weight	9 kg	Max. Load	130 kg
Lifting Column		Upholstery	Periline
- Height adjustment		Standard color	black
easyPro F & easyPro H	50 – 70 cm		
easyPro R	52 – 70 cm	Subject to technical modifications.	

easyPro

