



PRODUCT COMPARISON

performanceLine

UFSK
INTERNATIONAL

SYS



ophthalmic surgery



surgiTrend

Item No.: 2200 0012

surgiLine

Item No.: 2200 0008



chassis

zinc undercoated, powder-coated steel frame

chassis / dimensions

470 x 550 mm

470 x 600 mm

castor

5

4

brake system

central brake - 99% braking power on 4 castors

battery

strong & environmentally accumulators

mains operation

No

footswitch

height adjustment lifting column

top frame

zinc undercoated, powder-coated

seating systems

Item No.: 2200 0012

Comfort Sitz „Periline“

Item No.: 2200 0013

Item No.: 2200 0010

Dynamic Sitz „Stamskin Top“

Formsitz „Tundra skai“

Item No.: 2200 0008

height adjustment

Comfort seat 540 - 740 mm
Dynamic seat 600 - 800 mm

Form seat 560 - 760 mm
Comfort seat 540 - 740 mm

seat slide

No

armrests

stable, ergonomic shape, supporting area (350 x 120 mm)
height-adjustable, ball joint for full 360-degree range movement

height adjustment armrests

manual, asynchronous, non-sterile , 300 mm

**accessory armrest pair
for neurosurgery**

No

**armrest pair
with linear guide**

Yes (optional)
Item No.: 3351 1006

**armrest pair
with one-hand release**

Yes (optional)
Item No.: 3351 1004

footswitch plate adapter

Item No.: 3351 1011

Yes (optional)

Item No.: 3351 1014



AccutoMex