

<u>BRAND</u>	<u>MODEL</u>	<u>Merlin RA – Mounting</u>	<u>Merlin CLA – Mounting</u>	<u>ROLS[∞] / VROLS</u>
Alcon	LuxOR	Alcon Mounting Kit (11732)	Alcon Mounting Kit (11732)	manual (11313), auto (11363), VROLS-Z
Huvitz	HOM-700	Leica Mounting Kit (11409) & Vertical Spacing Kit (11425) (2 Spacers)	Leica Mounting Kit (11409) & Vertical Spacing Kit (11425) (2 Spacers)	manual (11315), auto (11364), VROLS-L
Inami	L-0980A	Inami Mounting Kit (11654)	Inami Mounting Kit (11654)	manual (11313), auto (11363), VROLS-Z
Leica	M220	DOES NOT MOUNT	DOES NOT MOUNT	VROLS-LX
Leica	M500	Leica Mounting Kit (11409)	Leica Mounting Kit (11409)	manual (11315), auto (11364), VROLS-L
Leica	M501	Leica Mounting Kit (11409)	Leica Mounting Kit (11409)	manual (11315), auto (11364), VROLS-L
Leica	M520	Leica Mounting Kit (11409); must remove counterweight plate	Leica Mounting Kit (11409); must remove counterweight plate	manual (11315), auto (11364), VROLS-L
Leica	M525	Leica Mounting Kit (11409); must remove counterweight plate	Leica Mounting Kit (11409); must remove counterweight plate	manual (11315), auto (11364), VROLS-L
Leica	M620	Leica Mounting Kit (11409)	Leica Mounting Kit (11409)	manual (11315), auto (11364), VROLS-L
Leica	M650	Leica Mounting Kit (11409)	Leica Mounting Kit (11409)	manual (11315), auto (11364), VROLS-L
Leica	M651	Leica Mounting Kit (11409)	Leica Mounting Kit (11409)	manual (11315), auto (11364), VROLS-L
Leica	M655	Leica Mounting Kit (11409)	Leica Mounting Kit (11409)	manual (11315), auto (11364), VROLS-L
Leica	M690	Leica 690 Mounting Kit for RA (11411)	Leica 690 Mounting Kit for CLA (11481)	manual (11315), auto (11364), VROLS-L
Leica	M691	Leica 690 Mounting Kit for RA (11411)	Leica 690 Mounting Kit for CLA (11481)	manual (11315), auto (11364), VROLS-L

<u>BRAND</u>	<u>MODEL</u>	<u>Merlin RA – Mounting</u>	<u>Merlin CLA – Mounting</u>	<u>ROLS[∞] / VROLS</u>
Leica	M695	Leica 690 Mounting Kit for RA (11411)	Leica 690 Mounting Kit for CLA (11481)	manual (11315), auto (11364), VROLS-L
Leica	M820	Leica Mounting Kit (11409) & Vertical Spacing Kit (11425)	Leica Mounting Kit (11409) & Vertical Spacing Kit (11425)	manual (11315), auto (11364), VROLS-L; mount above beam splitter
Leica	M822	Leica 690 Mounting Kit for RA (11411) & Vertical Spacing Kit (11425)	Leica 690 Mounting Kit for CLA (11481) & Vertical Spacing Kit (11425)	manual (11315), auto (11364), VROLS-L; mount above beam splitter
Leica	M840	Leica Mounting Kit (11409) & Vertical Spacing Kit (11425)	Leica Mounting Kit (11409) & Vertical Spacing Kit (11425)	manual (11315), auto (11364), VROLS-L; mount above beam splitter
Leica	M841	Leica Mounting Kit (11409) & Vertical Spacing Kit (11425)	Leica Mounting Kit (11409) & Vertical Spacing Kit (11425)	manual (11315), auto (11364), VROLS-L; mount above beam splitter
Leica	M844	Leica Mounting Kit (11409) & Vertical Spacing Kit (11425)	Leica Mounting Kit (11409) & Vertical Spacing Kit (11425)	manual (11315), auto (11364), VROLS-L; mount above beam splitter
Leica	Proveo 8	Leica Mounting Kit (11409) & Vertical Spacing Kit (11425)	Leica Mounting Kit (11409) & Vertical Spacing Kit (11425)	manual (11315), auto (11364), VROLS-L; mount above beam splitter
Möller-Wedel (Haag-Streit)	Allegra 90	DOES NOT MOUNT	DOES NOT MOUNT	manual (11313), auto (11363), VROLS-Z
Möller-Wedel (Haag-Streit)	Allegra 590	DOES NOT MOUNT	DOES NOT MOUNT	manual (11313), auto (11363), VROLS-Z
Möller-Wedel (Haag-Streit)	Allegra 900	DOES NOT MOUNT	DOES NOT MOUNT	manual (11313), auto (11363), VROLS-Z
Möller-Wedel (Haag-Streit)	CooperVision	Moeller Wedel Mounting Kit (11410)	Moeller Wedel Mounting Kit (11410)	manual (11313), auto (11363), VROLS-Z
Möller-Wedel (Haag-Streit)	Hi-R 900	Moeller Wedel Mounting Kit (11410)	Moeller Wedel Mounting Kit (11410)	manual (11313), auto (11363), VROLS-Z
Möller-Wedel (Haag-Streit)	Hi-R 900 SL	DOES NOT MOUNT	DOES NOT MOUNT	manual (11313), auto (11363), VROLS-Z
Möller-Wedel (Haag-Streit)	Hi-R 900 SLA	DOES NOT MOUNT	DOES NOT MOUNT	manual (11313), auto (11363), VROLS-Z

<u>BRAND</u>	<u>MODEL</u>	<u>Merlin RA – Mounting</u>	<u>Merlin CLA – Mounting</u>	<u>ROLS[∞] / VROLS</u>
Möller-Wedel (Haag-Streit)	EOS 900	DOES NOT MOUNT	DOES NOT MOUNT	manual (11313), auto (11363), VROLS-Z
Möller-Wedel (Haag-Streit)	Hi-R NEO 900	Moeller Wedel Mounting Kit (11410)	Moeller Wedel Mounting Kit (11410)	manual (11313), auto (11363), VROLS-Z
Möller-Wedel (Haag-Streit)	Ophthalmic 900	DOES NOT MOUNT	DOES NOT MOUNT	manual (11313), auto (11363), VROLS-Z
Möller-Wedel (Haag-Streit)	Spectra 900	DOES NOT MOUNT	DOES NOT MOUNT	manual (11313), auto (11363), VROLS-Z
Prescott	OMNI-FLEX	Moeller Wedel Mounting Kit (11410)	Moeller Wedel Mounting Kit (11410)	manual (11313), auto (11363), VROLS-Z
Storz	Urban	DOES NOT MOUNT	DOES NOT MOUNT	DOES NOT MOUNT
Takagi	OM-8	DOES NOT MOUNT	DOES NOT MOUNT	manual (11313), auto (11363), VROLS-Z
Takagi	OM-9	Takagi Mounting Kit (11412) and Takagi O06-25	Takagi Mounting Kit (11412) and Takagi O06-25	manual (11313), auto (11363), VROLS-Z
Takagi	OM-18	Takagi Mounting Kit (11412)	Takagi Mounting Kit (11412)	manual (11313), auto (11363), VROLS-Z
Topcon	OMS-90	Topcon OMS-90 Mouting Kit (11947)	Topcon OMS-90 Mouting Kit (11947)	manual (11313), auto (11363), VROLS-Z
Topcon	OM-600	DOES NOT MOUNT	DOES NOT MOUNT	manual (11313), auto (11363), VROLS-Z
Topcon	OMS-710	Topcon Mounting Kit (11413)	Topcon Mounting Kit (11413)	manual (11313), auto (11363), VROLS-Z
Topcon	OMS-800 Standard	Topcon Mounting Kit (11413)	Topcon Mounting Kit (11413)	manual (11313), auto (11363), VROLS-Z
Topcon	OMS-800 Pro	Topcon Mounting Kit (11413)	Topcon Mounting Kit (11413)	manual (11313), auto (11363), VROLS-Z

<u>BRAND</u>	<u>MODEL</u>	<u>Merlin RA – Mounting</u>	<u>Merlin CLA – Mounting</u>	<u>ROLS[∞] / VROLS</u>
Zeiss	Lumera	Zeiss Mounting Kit (11408)	Zeiss Mounting Kit (11408)	manual (11313), auto (11363), VROLS-Z; mount above beam splitter
Zeiss	Lumera T	Zeiss Mounting Kit (11408)	Zeiss Mounting Kit (11408)	manual (11313), auto (11363), VROLS-Z; mount above beam splitter
Zeiss	Lumera 300 Basic	DOES NOT MOUNT	DOES NOT MOUNT	manual (11313), auto (11363), VROLS-Z
Zeiss	Lumera 300 Standard	<i>field data needed</i>	<i>field data needed</i>	manual (11313), auto (11363), VROLS-Z
Zeiss	Lumera 300 Advanced	Zeiss Mounting Kit (11408)	Zeiss Mounting Kit (11408)	manual (11313), auto (11363), VROLS-Z
Zeiss	MDI	Zeiss Mounting Kit (11408)	Zeiss Mounting Kit (11408)	manual (11313), auto (11363), VROLS-Z
Zeiss	MDO	ZeissMD Mounting Kit for RA (11552)	ZeissMD Mounting Kit for CLA (11553)	manual (11313), auto (11363), VROLS-Z
Zeiss	MDU	ZeissMD Mounting Kit for RA (11552)	ZeissMD Mounting Kit for CLA (11553)	manual (11313), auto (11363), VROLS-Z
Zeiss	OPMI 1	Zeiss Mounting Kit (11408)	Zeiss Mounting Kit (11408)	manual (11313), auto (11363), VROLS-Z
Zeiss	OPMI 6	Zeiss Mounting Kit (11408)	Zeiss Mounting Kit (11408)	manual (11313), auto (11363), VROLS-Z
Zeiss	OPMI 6-CFC	Zeiss Mounting Kit (11408)	Zeiss Mounting Kit (11408)	manual (11313), auto (11363), VROLS-Z
Zeiss	OPMI 6S	Zeiss Mounting Kit (11408)	Zeiss Mounting Kit (11408)	manual (11313), auto (11363), VROLS-Z
Zeiss	OPMI 6-SFC	Zeiss Mounting Kit (11408)	Zeiss Mounting Kit (11408)	manual (11313), auto (11363), VROLS-Z
Zeiss	OPMI 6-SFR	Zeiss Mounting Kit (11408)	Zeiss Mounting Kit (11408)	manual (11313), auto (11363), VROLS-Z

BRAND	MODEL	Merlin RA – Mounting	Merlin CLA – Mounting	ROLS[∞] / VROLS
Zeiss	OPMI 6-CH	Zeiss Mounting Kit (11408)	Zeiss Mounting Kit (11408)	manual (11313), auto (11363), VROLS-Z
Zeiss	OPMI CS shell <i>(w/ various modules)</i>	<i>See module information below</i>	<i>See module information below</i>	manual (11313), auto (11363), VROLS-Z
Zeiss	Retrolux CS <i>on OPMI CS base</i>	Zeiss Mounting Kit (11408)	Zeiss Mounting Kit (11408)	<i>see OPMI CS shell information</i>
Zeiss	Retrolux 2 <i>on OPMI CS base</i>	DOES NOT MOUNT	DOES NOT MOUNT	<i>see OPMI CS shell information</i>
Zeiss	Retrolux 3 <i>on OPMI CS base</i>	DOES NOT MOUNT	DOES NOT MOUNT	<i>see OPMI CS shell information</i>
Zeiss	Retroskop <i>on OPMI CS base</i>	DOES NOT MOUNT	DOES NOT MOUNT	<i>see OPMI CS shell information</i>
Zeiss	Retroskop CS <i>on OPMI CS base</i>	<i>field data needed</i>	<i>field data needed</i>	<i>see OPMI CS shell information</i>
Zeiss	OPMI Pico	Zeiss Mounting Kit (11408); needs M4 x 12mm screws	Zeiss Mounting Kit (11408); needs M4 x 12mm screws	manual (11313), auto (11363), VROLS-Z
Zeiss	OPMI Vario	Zeiss Mounting Kit (11408)	Zeiss Mounting Kit (11408)	manual (11313), auto (11363), VROLS-Z
Zeiss	OPMI VISU 140	DOES NOT MOUNT	DOES NOT MOUNT	manual (11313), auto (11363), VROLS-Z
Zeiss	OPMI VISU 150	Zeiss Mounting Kit (11408)	Zeiss Mounting Kit (11408)	manual (11313), auto (11363), VROLS-Z
Zeiss	OPMI VISU 160	Zeiss Mounting Kit (11408)	Zeiss Mounting Kit (11408)	manual (11313), auto (11363), VROLS-Z
Zeiss	OPMI VISU 200	Zeiss Mounting Kit (11408)	Zeiss Mounting Kit (11408)	manual (11313), auto (11363), VROLS-Z
Zeiss	OPMI VISU 210	Zeiss Mounting Kit (11408)	Zeiss Mounting Kit (11408)	manual (11313), auto (11363), VROLS-Z

Last update: 30 March 2017