

# surgiLine

*The modern surgeon's stool.*

Distribuido por: Accutomex SA de CV  
www.accutomex.com  
ventas@accutomex.com  
Tel. +52 (81) 8347-2520 / 01-800-800-0507

Order No. 22000008

## ADVANTAGES



### Armrests

- Rugged, 360-degree movable, and height-adjustable
- Designed for personalized adjustment to the surgeon

### Seat System

- Seat system evenly distributes your weight as you make the slightest adjustment in position.
- Anatomically formed, ergonomic backrest reduces pressure on inter vertebral discs and provides ideal support for the lumbar spine.
- Seat system upholstery in standard color anthracite F6461197
- Outstanding sitting comfort

### Chassis

- 4 double casters
- Central brake acting on 4 wheels
- Power supply-independent battery box with integrated emergency stop button
- Foot keypad for moving the seat up & down
- Extremely easy to move
- Space-saving



## EQUIPMENT




Completely height and lengthwise adjustable arm-rests, 350 x 120 mm.



Optional: one-hand height release for direct and high-precision all-around adjustment.  
No. 33100000



Adjustable backrest.



Tiltable seat.



Foot keypad enables quick and easy electromotive height adjustment.



Anatomically formed, ergonomic backrest reduces pressure on intervertebral discs and provides ideal support for the lumbar spine.

Seat system evenly distributes your weight as you make the slightest adjustment in position.



Optional: adaptable foot switch plate.  
Order No. 33510014

## SPECIFICATIONS

### Lifting Column

- Drive 3,000 N IPX4
- Lifting height 200 mm

### Chassis

- w x l 470 x 535 mm
- Wheels ø 100 mm
- Central brake acting on all four wheels
- Emergency stop button

Weight 45 kg

Max Load 150 kg

### Battery unit

- Battery box with 2 batteries (lead-gel), pollution free,

maintenance-free, cadmium-free

2 x 12 V 7.2 Ah

- Universal battery charger:

Input: 100-240 V

Output: 24 V

### Seat/backrest

- Min. sitting height 535 mm

- Max. sitting height 735 mm

- Seat adjusting +6°/-9°

- Back adjusting +4°/-18°

Upholstery colors Tundra skai

Standard color anthracite

F6461197

Subject to technical modifications.

